 **Manufacturer:**

REGEDENT AG  
Zollikerstrasse 144  
CH-8008 Zürich  
Switzerland  
www.regedent.com

 **REGEDENT**  
smart regeneration

CE

# PERISOLV

INSTRUCTIONS FOR USE

**EN** ENGLISH

## DEVICE DESCRIPTION

Oral Debridement Gel

One PERISOLV unit of use consists of a sealed and labelled pouch containing:

- one syringe filled with an alkaline 0.95 % sodium hypochlorite solution, closed with a cap
- one syringe filled with a slightly viscous alkaline gel containing glutamic acid, leucine, lysine, carboxymethyl cellulose and titanium dioxide, closed with a cap
- one connector to connect the two syringes and mix the contents
- one non-invasive blunt applicator tip for application of the mixture

The mixing of the two components just before application results in an opaque gel containing 0.43 % sodium hypochlorite with a slight odour of chlorine. The gel is applied to the treatment site using the applicator tip. After removal of the gel, debridement is performed to remove biofilm from the treatment site. Five (5) units of use are packed together in one labelled cardboard shelf box.

## INDICATION FOR USE

PERISOLV is indicated for patients diagnosed with periodontitis, mucositis or peri-implantitis.

## INTENDED USE

PERISOLV is intended to be used non-surgically as a supplementary treatment, applied before debridement, to precondition the treatment site to enable oral tissue regeneration. PERISOLV is a gel intended for softening and degrading biofilm in periodontal and peri-implant pockets.

## INTENDED PATIENT POPULATION

PERISOLV can be used for all patients considered fit for treatment by dental professional performing the treatment. Special considerations shall be given to:

- Paediatric patients
- Nursing mothers
- Pregnant women
- Patients allergic to PERISOLV ingredients

EN

DE

ES

FR

IT

NL

DK

SE

NO

FI

HU

PL

RU

RO

## **INSTRUCTION FOR USE**

### **Pre-treatment**

1. Plan to treat one quadrant at a time.
2. Keep the product at room temperature a short time before use. Take the syringes out of the aluminium pouch together with the connector and the blunt applicator tip. Hold the syringe with the opening upwards and remove the caps.
3. Apply the connector and mount the syringes together.
4. Mix the two components by pressing alternately on the end of each syringe 10–15 times until the solution is homogenous.
5. Press all the liquid into the transparent syringe. Remove the connector. Attach the blunt applicator tip to the transparent syringe.

### **Treatment**

6. Immediately after mixing, apply PERISOLV to the site using the blunt applicator tip. When treating a pocket, fill with PERISOLV up to the gum line starting at the most apical aspect.
7. Leave PERISOLV for at least 30 seconds on-site, but no longer than 5 minutes, before rinsing with saline solution.
8. Perform debridement according to routine procedure.
9. If necessary, repeat steps 6 to 8.

## **INTENDED USERS**

Dental practitioners, such as professional dentists, and dental hygienists, experienced in treating periodontitis and/or peri-implantitis, as well as mucositis. The product is not intended to be used by a layperson or any other person outside the intended users.

## **INTENDED USE ENVIRONMENT**

The product is intended to be used in an adequately sterile field, such as an operating or treatment room at a dental office, or in a dental clinic.

## **PRECAUTIONS**

- PERISOLV must not be used in combination with chlorhexidine, hydrogen peroxide or strong acids.
- In case of contact of PERISOLV with the eyes, rinse carefully with plenty of water for several minutes.
- 0.95 % sodium hypochlorite solution and PERISOLV can stain clothing.
- Do not use PERISOLV after the expiry date indicated on the label or when the aluminium pouch is damaged.

## **HANDLING, STORAGE AND DISPOSAL**

For long term storage keep the product refrigerated at 2–8 °C / 36–46 °F in its original packaging. Do not freeze. Short term exposure of cumulatively 72 hours to temperatures of maximum 25 °C / 77 °F will not affect the product performance.

Disposal should be handled in accordance with local regulations and instructions.

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